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### Preparation of the patient

#### 1. Cleaning the tooth surface

Clean the tooth surface thoroughly by the usual method. Pay attention to the oil-free components and complete removal of substances by rinsing and subsequent drying of the tooth with compressed air.

#### 2. Conditioning of the tooth surface

Condition the entire tooth surface and leave the acid gel in place for the time specified in the manufacturer's instructions. Remove the gel by aspiration and rinsing. When using a self-conditioning primer, apply the primer to the tooth surface and let it act.

#### 3. Isolation of the oral cavity

The brackets can be fixed without any additional preparations. Ensure that the bonding area is dry. Saliva and water can prevent successful bonding of the brackets.


#### 4. Applying the primer (bonding agent)

Now apply a thin layer of primer (bonding agent) to the tooth surface and gently dry with compressed air. If using a self-conditioning primer, omit this step.


## Indirect bonding with INDIVIDUA®

### Placing the brackets

#### 1. Applying adhesive to the bracket base

Apply a thin layer of primer to the bracket base, then apply the adhesive, in correct portions. It is recommended to spread the adhesive towards the apical end of the bracket base, where excess can be removed. Ensure the the correct amount of adhesive.

#### 2. Smoothing the adhesive on the bracket base

Use an instrument wetted with primer to spread the adhesive very thinly over the bracket base, ensuring that the whole surface is evenly covered.

#### 3. Placement of trays

Ensure that you position and place the trays exactly on the correct teeth. Avoid contact between the prepared bracket bases and the adjacent teeth. Push the tray gently in the apical direction until it is precisely in position. Press against the tray base to check that it is not loose.

#### 4. Removal of excess adhesive

Excess adhesive can escape at the apical edge of the bracket base. Run an instrument along the apical edge of the bracket base and remove the excessive adhesive. Carefully check this area once again before it cures.

#### 5. UV curing

Cure the adhesive with a high-power LED lamp. Always check the required curing time for the system used.

#### 6. Tray removal

Remove the trays using an instrument placed apically at the ends of the trays at the bottom edge. Pull the tray gently towards occlusal, alternating between the ends of the tray. Secure the tray so that it cannot be swallowed.

#### 7. End result

Now all the brackets are fixed and you can proceed with your orthodontic treatment.